Cat. No. BME100768



## **PRODUCT INFORMATION**

**Uniprot ID** P18509

LY3451838;LY-3451838;LY 3451838 **Common Name** 

Conjugate Unconjugated

**PACAP Synonyms** N/A **Applications** 

Recommended

N/A **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution for specific instructions of reconstitution.

**Host Species** N/A IgG4 IgG type Reactivity Human **Target** ADCYAP1

**Description** Anti-ADCYAP1(vipalanebart biosimilar IgG4) mAb

**Delivery** 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store

> > Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only

