

PRODUCT INFORMATION

Uniprot ID P05067

Common Name LY-3372993;LY 3372993;LY3372993;N3pG 4;N3pG4;N3PG-IV

Conjugate Unconjugated

AAA;AD1;PN2;ABPP;APPI;CVAP;ABETA;PN-II;preA4;CTFgamma;alpha-sAPP **Synonyms**

Applications Recommended N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species IgG type lgG1 Reactivity Human APP **Target**

Anti-APP(remternetug biosimilar IgG1) mAb Description

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. Storage & Shipping

Background

Email: info@dimabio.com Website: www.dimabio.com

Research use only Usage

