

## **PRODUCT INFORMATION**

**Uniprot ID** P30530

BGB 149;BGB-149;BGB149 **Common Name** 

Conjugate Unconjugated

AXL;ARK;UFO;AXL3;JTK11;Tyro7 **Synonyms** 

**Applications** N/A

Recommended

N/A **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

**Host Species** N/A lgG1 IgG type Reactivity Human **Target AXL** 

**Description** Anti-AXL(tilvestamab biosimilar IgG1) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only





