Cat. No. BME100744



## **PRODUCT INFORMATION**

**Uniprot ID** C7G0W1

ZTS-00008183;ZTS00008183;ZTS **Common Name** 

00008183;ZTS-8183;ZTS8183;ZTS 8183

Conjugate Unconjugated

Interleukin-31;IL-31 **Synonyms** 

**Applications** N/A

Recommended Dilutions

N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

**Host Species** IgG2 IgG type Reactivity Canine

Canine IL31 **Target** 

Anti-Canine IL31(riltovetbart biosimilar IgG2) Description

mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

temperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only

