

## **PRODUCT INFORMATION**

**Uniprot ID** P51685

**Common Name** BAY 3375968;BAY-3375968;BAY3375968

Conjugate Unconjugated

**Synonyms** CY6;TER1;CCR-8;CKRL1;CDw198;CMKBR8;GPRCY6;CMKBRL2;CC-CKR-8

**Applications** Recommended N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

**Host Species** N/A IgG type lgG1 Reactivity Human CCR8 **Target** 

Description Anti-CCR8(lanerkitug biosimilar [IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

Email: info@dimabio.com Website: www.dimabio.com

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background** 

Research use only Usage

