Cat. No. BME100597



PRODUCT INFORMATION

Uniprot ID P11836

Common Name Arzerra; GSK-1841157; GSK 1841157; GSK

1841157;Kesimpta;OMB-157;OMB157;OMB 157

Conjugate Unconjugated

Synonyms B1;S7;Bp35;CD20;FMC7;CVID5;LEU-16;MS4A1

Applications N/A
Recommended N/A

Dilutions IN//

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A
IgG type IgG1
Reactivity Human
Target CD20

Description Anti-CD20(ofatumumab biosimilar [IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

Storage & Shipping intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

