

## **PRODUCT INFORMATION**

**Uniprot ID** P28907

GEN-3014; GEN3014; GEN 3014 **Common Name** 

Conjugate Unconjugated

CD38;ADPRC1;cADPR1;ADPRC 1;T10;cADPr **Synonyms** 

hydrolase 1

**Applications** N/A

Recommended

N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

**Host Species** lgG1 IgG type Reactivity Human **CD38 Target** 

**Description** Anti-CD38(erzotabart biosimilar IgG1) mAb

4-6 weeks **Delivery** 

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage** 

