

## **PRODUCT INFORMATION**

**Uniprot ID** P28907

MOR 202; MOR202; MOR-202; MOR-03087; MOR

03087;MOR03087;MOR **Common Name** 

N/A

3087; MOR-3087; MOR3087; TJ-202; TJ202; TJ 202

Conjugate Unconjugated

CD38;ADPRC1;cADPR1;ADPRC 1;T10;cADPr **Synonyms** 

hydrolase 1

**Applications** N/A

Recommended

**Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

**Host Species** N/A lgG1 IgG type Reactivity Human **Target CD38** 

Description Anti-CD38(felzartamab biosimilar ☐ lgG1) mAb

4-6 weeks **Delivery** 

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

temperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage** 

