

PRODUCT INFORMATION

Uniprot ID P28907

Common Name hu38SB19;SAR-650984;SAR650984;Sarclisa

Conjugate Unconjugated

CD38;ADPRC1;cADPR1;ADPRC 1;T10;cADPr **Synonyms**

hydrolase 1

Applications N/A

Recommended

N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG1 IgG type Reactivity Human **CD38 Target**

Description Anti-CD38(isatuximab biosimilar ☐ IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

Email: info@dimabio.com Website: www.dimabio.com

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

