Cat. No. BME100499



## **PRODUCT INFORMATION**

**Uniprot ID** P49961

ES 002;ES002;ES-002;ES 002023;ES-002023;ES **Common Name** 

002023

Unconjugated Conjugate

CD39:SPG64:ATP-DPH:ATPDase:NTPDase-1:ENTPD1 **Synonyms** 

**Applications** N/A

Recommended

N/A **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % Formulation & - 8% trehalose is added as protectants before

Reconstitution lyophilization. Please see Certificate of Analysis for

specific instructions of reconstitution.

**Host Species** N/A IgG type IgG1 Reactivity Human **Target CD39** 

**Description** Anti-CD39(eurestobart biosimilar IgG1) mAb

**Delivery** 4-6 weeks

Email: info@dimabio.com Website: www.dimabio.com

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid Storage & Shipping

repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans or

animals.

Usage Research use only

