

PRODUCT INFORMATION

Uniprot ID P07766

SPV-T3a;SPVT3a **Common Name** Conjugate Unconjugated

Synonyms T3E;TCRE;IMD18;CD3epsilon

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG2b IgG type Reactivity Human **Target** CD3E

Description Anti-CD3E(dafsolimab biosimilar IgG2b) mAb

Delivery 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store

Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only



