

## **PRODUCT INFORMATION**

**Uniprot ID** P29965

**Common Name** AT-1501;AT1501;AT 1501

Conjugate Unconjugated

IGM;IMD3;TRAP;gp39;CD154;CD40L;HIGM1;T-BAM;TNFSF5;hCD40L **Synonyms** 

N/A **Applications** Recommended Dilutions N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

**Host Species** IgG type lgG1 Reactivity Human **Target** CD40 Ligand

Anti-CD40 Ligand(tegoprubart biosimilar ☐ IgG1) mAb **Description** 

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. Storage & Shipping

Email: info@dimabio.com Website: www.dimabio.com

**Background** 

Usage Research use only

