Cat. No. BME100573



PRODUCT INFORMATION

Uniprot ID P25942

ADC 1013;ADC-1013;ADC1013;JNJ 7107;JNJ7107;JNJ-7107;JNJ-64457107;JNJ64457107;JNJ64457107;JNJ **Common Name**

Conjugate Unconjugated

CD40;p50;Bp50;CDW40;TNFRSF5 **Synonyms**

Applications N/A Recommended

Dilutions

N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Formulation & Reconstitution

specific instructions of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human **CD40 Target**

Description Anti-CD40(mitazalimab biosimilar IgG1) mAb

Delivery

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized

form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid

repeated freezing and thawing). Lyophilized proteins

are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic **Background**

or diagnostic procedures for humans or animals.

Email: info@dimabio.com Website: www.dimabio.com

Research use only Usage

