

PRODUCT INFORMATION

Uniprot ID P25942

ABBV 323;ABBV-323;ABBV 323;PR **Common Name** 1629977;PR-1629977;PR1629977

Conjugate Unconjugated

CD40;p50;Bp50;CDW40;TNFRSF5 **Synonyms**

Applications N/A Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG1 IgG type Reactivity Human **CD40 Target**

Description Anti-CD40(ravagalimab biosimilar IgG1) mAb

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**



Email: info@dimabio.com Website: www.dimabio.com

