

PRODUCT INFORMATION

|                              |   |
|------------------------------|---|
| Uniprot ID                   | P25942  |
| Common Name                  | APX 005;APX005;APX 005;APX5;APX-5;APX 5;EPI-0050;EPI0050;EPI 0050;EPI 50;EPI-50;EPI50   |
| Conjugate                    | Unconjugated  |
| Synonyms                     | CD40;p50;Bp50;CDW40;TNFRSF5   |
| Applications                 | N/A   |
| Recommended Dilutions        | N/A   |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.  |
| Host Species                 | N/A   |
| IgG type                     | IgG1  |
| Reactivity                   | Human   |
| Target                       | CD40  |
| Description                  | Anti-CD40(sotigalimab biosimilar IgG1) mAb  |
| Delivery                     | 4-6 weeks   |
| Storage&Shipping             | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background                   | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.   |
| Usage                        | Research use only   |

