

PRODUCT INFORMATION

Uniprot ID P25942

APX 005;APX005;APX 005;APX5;APX-5;APX **Common Name**

5;EPI-0050;EPI0050;EPI 0050;EPI 50;EPI-50;EPI50

Conjugate Unconjugated

CD40;p50;Bp50;CDW40;TNFRSF5 **Synonyms**

Applications N/A Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species IgG1 IgG type Reactivity Human **CD40 Target**

Description Anti-CD40(sotigalimab biosimilar IgG1) mAb

Delivery

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

