Cat. No. BME100464



PRODUCT INFORMATION

Uniprot ID P32970

Common Name

ARGX-110;ARGX110;JNJ 74494550;JNJ74494550;JNJ-74494550;JNJ-4550;JNJ4550

Conjugate Unconjugated

Synonyms CD27L;LPFS3;CD27-L;CD27LG;TNFSF7;TNLG8A

Applications N/A Recommended

N/A

Formulation & Reconstitution

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % -8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for

specific instructions of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human **Target CD70**

Anti-CD70(cusatuzumab biosimilar IgG1) mAb **Description**

Delivery 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid

> > Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping repeated freezing and thawing). Lyophilized proteins

are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic **Background**

or diagnostic procedures for humans or animals.

Usage Research use only

