Cat. No. BME100627



PRODUCT INFORMATION

Uniprot ID P06731

M9140;M-9140;M 9140 **Common Name**

Conjugate Unconjugated

CEACAM-5;CEA;CD66e;Meconium antigen 100 **Synonyms**

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before Formulation & Reconstitution lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A lgG1 IgG type Reactivity Human CEACAM5 **Target**

Anti-CEACAM5(precemtabart biosimilar IgG1) **Description**

mAb

Delivery 4-6 weeks

> Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

therapeutic or diagnostic procedures for humans **Background**

or animals.

Research use only **Usage**

