

PRODUCT INFORMATION

Uniprot ID P56856

Common Name LM-102;LM 102;LM102

Conjugate Unconjugated **Synonyms** SFTA5;SFTPI

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A lgG1 IgG type Reactivity Human **Target** CLDN18

Description Anti-CLDN18(tecotabart biosimilar ☐ IgG1) mAb

Delivery 4-6 weeks

Storage & Shipping

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com


témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only



/+86-400-006-0995(China)

