

PRODUCT INFORMATION

|                              |                                                                                                                                                                                                                                                     |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Uniprot ID                   | P16410                                                                                                                                                                                                                                              |
| Common Name                  | BMS-986249;BMS 986249;BMS986249                                                                                                                                                                                                                     |
| Conjugate                    | Unconjugated                                                                                                                                                                                                                                        |
| Synonyms                     | CD;GSE;GRD4;ALPS5;CD152;CTLA-4;IDDM12;CELIAC3                                                                                                                                                                                                       |
| Applications                 | N/A                                                                                                                                                                                                                                                 |
| Recommended Dilutions        | N/A                                                                                                                                                                                                                                                 |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.                                                    |
| Host Species                 | N/A                                                                                                                                                                                                                                                 |
| IgG type                     | IgG1                                                                                                                                                                                                                                                |
| Reactivity                   | Human                                                                                                                                                                                                                                               |
| Target                       | CTLA-4                                                                                                                                                                                                                                              |
| Description                  | Anti-CTLA-4(sovipostobart biosimilar-IgG1) mAb                                                                                                                                                                                                      |
| Delivery                     | 3-4 weeks                                                                                                                                                                                                                                           |
| Storage & Shipping           | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background                   | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.                                                                                                                                               |
| Usage                        | Research use only                                                                                                                                                                                                                                   |

