Cat. No. BME100612



PRODUCT INFORMATION

Uniprot ID P25101

GMA-301;GMA301;GMA 301;Getagozumab **Common Name**

Conjugate Unconjugated

ETA; ET-A; ETAR; ETRA; MFDA; ETA-R; hET-AR **Synonyms**

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG4 IgG type Reactivity Human **Target EDNRA**

Description Anti-EDNRA(patecibart biosimilar IgG4) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

