

PRODUCT INFORMATION

Uniprot ID P00533

DS-992;DS992;S-95026;S95206;Sym-004 (Combination with Modotuximab);Sym004 (Combination with Modotuximab) **Common Name**

Conjugate Unconjugated

Synonyms ERBB;ERRP;HER1;mENA;ERBB1;NNCIS;PIG61;NISBD2

Applications N/A Recommended **Dilutions**

N/A

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Host Species N/A IgG type lgG1 Reactivity Human **Target EGFR**

Anti-EGFR(futuximab biosimilar ☐ IgG1) mAb Description

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawning). Lyophilized proteins Storage & Shipping

are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic **Background**

or diagnostic procedures for humans or animals.

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only

