Cat. No. BME100702



## **PRODUCT INFORMATION**

**Uniprot ID** P00533

**Common Name** GT-MAB 5.2-GEX Conjugate Unconjugated

ERBB;ERRP;HER1;mENA;ERBB1;NNCIS;PIG61;NISBD2 **Synonyms** 

**Applications** N/A Recommended N/A

**Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Formulation & Reconstitution

specific instructions of reconstitution.

**Host Species** IgG type lgG1 Reactivity Human **Target EGFR** 

Anti-EGFR(tomuzotuximab biosimilar IgG1) mAb **Description** 

**Delivery** 4-6 weeks

Store at -20  $^{\circ}$ C to -80  $^{\circ}$ C for 12 months in lyophilized form. After reconstitution, if not intended for use

within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

Research grade biosimilar. Not for use in therapeutic **Background** 

or diagnostic procedures for humans or animals.

Email: info@dimabio.com Website: www.dimabio.com

**Usage** Research use only

