PRODUCT INFORMATION

Clone ID DM147 EPCAM Target

EPCAM;TACSTD1;TROP1;CD326;DIAR5;EGP2;EGP314;EGP40;ESA;GA733-2;HNPCC8;HNPCC-8;KS1:4;KSA;M4S1;MIC18;MK1 Synonyms

Host Species Rabbit

Description Anti-EPCAM antibody(DM147): Rabbit mAb

Delivery In Stock P16422 Uniprot ID IgG type Rabbit IgG Clonality Monoclonal Reactivity Human Applications ELISA; Flow Cyt

Recommended Dilutions ELISA 1:5000-10000: Flow Cvt 1:100

Purification Purified from cell culture supernatant by affinity chromatography

Formulation & Reconstitution Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage & Shipping

Steve et-20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. This gene encodes a carcinoma-associated antigen and is a member of a family that includes at least two type I membrane proteins. This antigen is expressed on most normal epithelial cells and gastrointestinal carcinomas and functions as a homotypic calcium-independent cell adhesion molecule. The antigen is being used as a target for immunotherapy treatment of human carcinomas. Mutations in this gene result in congenital tufting enteropathy. Background

Usage Research use only Conjugate Unconjugated

All DIMA recombinant antibodies are genuinely generated by DIMA Biotech. They are all under patent application. Any protein sequencing or reverse engineering attempt is prohibited. We are actively scrutinizing all patent application to ensure no IP infringement. **DIMA Disclaimer**




Figure 1. ELISA plate pre-coated by 1 μ g/ml (100 μ l/well) Human EPCAM protein, His tagged protein PME100068 can bind Rabbit anti-EPCAM monoclonal antibody (clone: DM147) in a linear range of 5-100 ng/ml.

Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)

Email: info@dimabio.com Website: www.dimabio.com



Cat. No. DME100147




Figure 2. EPCAM protein is highly expressed on the surface of HEK293 cell membrane. Flow cytometry analysis with Anti-EPCAM (DM147) on HEK293 cells transfected with human EPCAM (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram), and Isotype antibody on HEK293 transfected with irrelevant protein (Orange histogram).




Figure 3. Flow cytometry analysis of antigen binding of Rabbit anti-human EPCAM mAb(DME100147). (A) DME100147 does not bind to Jurkat cells that do not express EPCAM. (B) A clear peak shift of DME100147 was seen compared to the control when incubated with EPCAM -expressing HT55 cells, indicating strong binding of DME100147 to EPCAM. Antibodies were incubated at 10 μ g/mL.

Email: info@dimabio.com Website: www.dimabio.com



Address: Wuhan institute of Biotechnology B7, Biolake No.666 Gaoxin Road, Wuhan, Hubei, China Telephone: +1 2409940618(USA) /+86-18062749453(China) /+86-400-006-0995(China)