

PRODUCT INFORMATION

Uniprot ID	O14944
Common Name	LY3016859, LY-3016859
Conjugate	Unconjugated
Synonyms	EPR
Applications	Flow Cyt
Recommended Dilutions	Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Humanized
IgG type	Human IgG4(F234A,L235A) - Kappa
Reactivity	Human
Target	EREG
Description	Anti-EREG(fepixnebart biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only






Figure 1. Flow cytometry analysis with 7.5 µg/mL Anti-EREG(fepixnebart biosimilar) mAb (BME100511) on HEK293 cells transfected with Human EREG protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

