

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | FPA-144, FPA144, FPA114-A |
| Conjugate | Unconjugated |
| Synonyms | BEK, KGFR, KSAM |
| Applications | ELISA |
| Recommended Dilutions | ELISA 1:5000-10000 |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Humanized |
| IgG type | Human IgG1(K214R) - kappa |
| Reactivity | Human |
| Target | FGFR2 |
| Uniprot ID | P21802 |
| Description | Anti-FGFR2(bemarituzumab biosimilar) mAb |
| Delivery | In Stock |
| Storage&Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |




Anti-FGFR2(bemarituzumab biosimilar) mAb ELISA0.2 μ g of Human FGFR2, His tagged protein per well

Figure 1. ELISA plate pre-coated by 2 μ g/mL (100 μ L/well) Human FGFR2 Protein, His Tag (PME101408) can bind Anti-FGFR2(bemarituzumab biosimilar) mAb (BME100260) in a linear range of 0.13-16 ng/mL.

DIMABIO CONFIDENTIAL

