

## PRODUCT INFORMATION

**Uniprot ID** P02751

**Common Name** PYX-201; PYX201; PYX 201

Conjugate Unconjugated

**Synonyms** FN;CIG;FNZ;MSF;ED-B;FINC;GFND;LETS;GFND2;SMDCF

**Applications** N/A Recommended N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

**Host Species** lgG1 IgG type Reactivity Human FN1 **Target** 

Anti-FN1(micvotabart biosimilar IgG1) mAb **Description** 

**Delivery** 4-6 weeks

Storage & Shipping

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins

are shipped at ambient temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background** 

Email: info@dimabio.com Website: www.dimabio.com

Usage Research use only

