

PRODUCT INFORMATION

Uniprot ID Q99988

CTL-002;CTL002;CTL 002 **Common Name**

Conjugate Unconjugated

Synonyms HG;PDF;MIC1;PLAB;MIC-1;NAG-1;PTGFB;GDF-15

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG4 IgG type Reactivity Human **Target** GDF15

Description Anti-GDF15(visugromab biosimilar ☐ lgG4) mAb

Delivery 4-6 weeks

Storage & Shipping

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

> Email: info@dimabio.com Website: www.dimabio.com

Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

