

PRODUCT INFORMATION

Uniprot ID 014793

REGN-1033; REGN1033; REGN **Common Name**

1033;SAR-391786;SAR 391786;SAR391786

Conjugate Unconjugated

GDF8;MSLHP;MSTN **Synonyms**

Applications N/A

Recommended **Dilutions**

N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG4 IgG type Reactivity Human GDF8 **Target**

Description Anti-GDF8(trevogrumab biosimilar IgG4) mAb

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

