


PRODUCT INFORMATION

Clone ID	OTI2A4 (formerly 2A4)
Target	GST tag
Synonyms	ANTI-GST; ANTIBODY FOR GST; ANTI GST; GST AB; GST ANTIBODY
Host Species	Mouse
Description	Mouse monoclonal eGFP antibody, clone OTI5A2
Delivery	1-2week
Uniprot ID	N/A
IgG type	IgG1
Clonality	Monoclonal
Reactivity	N/A
Applications	WB
Recommended Dilutions	WB: 1:4000
Purification	N/A
Formulation & Reconstitution	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage&Shipping	Store at -20°C.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated



Anti-GST Antibody (2A4) Mouse mAb

Cat. No. DME900010



recombinant proteins, and one GST-DDK bacterial lysate produced in E.coli. at 1:5000 dilution

DIMABIO CONFIDENTIAL

Address: Wuhan institute of
Biotechnology B7, Biolake No.666
Gaoxin Road, Wuhan, Hubei, China
Telephone: +1 2409940618(USA)
/+86-18062749453(China)
/+86-400-006-0995(China)

Email: info@dimabio.com
Website: www.dimabio.com

