

PRODUCT INFORMATION

Clone ID	CB051
Target	HA tag
Synonyms	ANTI-HA; Anti-HA Antibody; Anti HA; Anti HA Antibody; ANTI INFLUENZA
Host Species	Mouse
Description	Carrier-free (BSA/glycerol-free) HA mouse monoclonal antibody, clone CB051
Delivery	3-4week
Uniprot ID	N/A
IgG type	IgG1
Clonality	Monoclonal
Reactivity	N/A
Applications	WB
Recommended Dilutions	WB: 1:2000
Purification	N/A
Formulation & Reconstitution	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Storage&Shipping	Store at -20°C.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL






Figure 1. HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or HA tagged LGALS3 cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HA.




Figure 2. Western Blot comparison of anti-HA Tag antibodies from DME900039, BioLegend (mms-101p), and Abcam (ab18181). 0.04ul (Lane 1, 3, 5) and 0.2ul (Lane 2, 4, 6) of HA-PARD6B lysate (recombinant proteins with HA tag at the N-terminus), blotted with 1:1000 dilution of TA180128 (lane 1-2), mms-101p (Lane 3-4), and ab18181 (Lane 5-6).

