

PRODUCT INFORMATION

Clone ID	CB051
Target	HA tag
Synonyms	ANTI-HA; Anti-HA Antibody; Anti HA; Anti HA Antibody; ANTI INFLUENZA
Host Species	Mouse
Description	Carrier-free (BSA/glycerol-free) HA mouse monoclonal antibody, clone CB051
Delivery	3-4week
Uniprot ID	N/A
IgG type	IgG1
Clonality	Monoclonal
Reactivity	N/A
Applications	WB
Recommended Dilutions	WB: 1:2000
Purification	N/A
Formulation & Reconstitution	Lyophilized powder (original buffer 1X PBS, pH 7.3, 8% trehalose)
Storage&Shipping	Store at -20°C.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated

DIMABIO CONFIDENTIAL

