Cat. No. BME100271

PRODUCT INFORMATION

Common Name Hen egg Lysozyme Conjugate Unconjugated

Synonyms HEL **Applications** N/A Recommended N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution for specific instructions of reconstitution.

Host Species Chimeric

IgG type Human IgG1(L234A,L235A) - Kappa

Reactivity Human **Target** HEL **Uniprot ID** N/A

Anti-HEL Human IgG1(L234A/L235A)-Kappa Isotype control mAb **Description**

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

Email: info@dimabio.com Website: www.dimabio.com

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage**

