

## **PRODUCT INFORMATION**

**Uniprot ID** P04626

PF-06741086;PF06741086;PF-6741086;PF6741086 **Common Name** 

Conjugate Unconjugated

NEU; NGL; HER2; TKR1; CD340; HER-2; VSCN2; MLN **Synonyms** 19;MLN-19;c-ERB2;c-ERB-2;HER-2/neu;p185(erbB2)

**Applications** 

Recommended

N/A **Dilutions** 

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for Formulation & Reconstitution

specific instructions of reconstitution.

**Host Species** N/A lgG1 IgG type Reactivity Human Her2 **Target** 

**Description** Anti-Her2(marstacimab biosimilar IgG1) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Storage & Shipping

Lyophilized proteins are shipped at ambient

temperature.

**Background** 

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans

Email: info@dimabio.com Website: www.dimabio.com

**Usage** Research use only

