

PRODUCT INFORMATION

Clone ID	1B9
Target	His Tag
Synonyms	his tag
Host Species	Rabbit
Description	Anti His (1B9) Chimeric mAb-Gallus IgY
Delivery	In Stock
Uniprot ID	N/A
IgG type	Rabbit/Gallus IgY
Clonality	Monoclonal
Reactivity	N/A
Applications	ELISA
Recommended Dilutions	ELISA 1/5000-10000
Purification	Purified from cell culture supernatant by affinity chromatography
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated






Figure 1. ELISA plate pre-coated by 1 µg/ml (100 µl/well) Mouse CDH17, His tagged protein can bind Rabbit anti-His monoclonal antibody(clone: 1B9) in a linear range of 5-100 ng/ml.

