

PRODUCT INFORMATION

Clone ID	OTI2B5 (formerly 2B5)
Target	His Tag
Synonyms	ANTI HIS, ANTI HISTIDINE, ANTIHIS
Host Species	Mouse
Description	Clone OTI2B5, Anti-His mouse monoclonal antibody
Delivery	1-2week
Uniprot ID	N/A
IgG type	IgG1
Clonality	Monoclonal
Reactivity	His tag
Applications	Flow CytWB
Recommended Dilutions	WB 1:4000
Purification	N/A
Formulation & Reconstitution	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Storage & Shipping	Store at -20°C.
Background	N/A
Usage	Research use only
Conjugate	Unconjugated

Figure 1. Anti-His Immunoblot analysis of six different HEK293T over-expression lysates containing different His tagged recombinant proteins, and two his tagged bacterial lysate produced in E.coli. at 1:4000 dilution.

