

PRODUCT INFORMATION

Uniprot ID Q9Y6W8

KY-1044;KY1044;SAR-445256 **Common Name**

Conjugate Unconjugated

AILIM;CD278;CVID1 **Synonyms**

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A lgG1 IgG type Reactivity Human **Target ICOS**

Description Anti-ICOS(alomfilimab biosimilar IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only



