

PRODUCT INFORMATION

| | |
|---|---|
| Common Name | AK101,AK-101 |
| Conjugate | Unconjugated |
| Synonyms | IL-12B;CLMF p40;IL-12 subunit p40;NKSF2 |
| Applications | ELISA |
| Recommended Dilutions | ELISA 1:5000-10000 |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Host Species | Humanized |
| IgG type | Human IgG1 - kappa |
| Reactivity | Human |
| Target | IL12B |
| Uniprot ID | P29460 |
| Description | Anti-IL12B(ebdarokimab biosimilar) mAb |
| Delivery | In Stock |
| Storage&Shipping | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature. |
| Background | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. |
| Usage | Research use only |




Anti-IL12B(ebdarokimab biosimilar) mAb ELISA0.2 μ g of Human IL12B, His tagged protein per well

Figure 1. ELISA plate pre-coated by 2 μ g/mL (100 μ L/well) Human IL12B Protein, His Tag (PME101133) can bind Anti-IL12B(ebdarokimab biosimilar) mAb (BME100223) in a linear range of 0.13-16 ng/mL.

