

PRODUCT INFORMATION

Common Name	MABp1
Conjugate	Unconjugated
Synonyms	IL-1 alpha;IL-1A;IL1;IL1-ALPHA;IL1F1
Applications	ELISA; Flow Cyt
Recommended Dilutions	ELISA 1:5000-10000; Flow Cyt 1:100
Formulation & Reconstitution	Lyophilized from sterile PBS, pH 7.4. Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Host Species	Humanized
IgG type	Human IgG1 - kappa
Reactivity	Human
Target	IL1A
Uniprot ID	P01583
Description	Anti-IL1A(bermekimab biosimilar) mAb
Delivery	In Stock
Storage&Shipping	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Background	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Usage	Research use only



Anti-IL1A (bermekimab biosimilar) mAb ELISA

0.2 μ g of Human IL1A, hFc tagged protein per well




Figure 1. ELISA plate pre-coated by 2 μ g/mL (100 μ L/well) Human IL1A Protein, hFc Tag (PME101053) can bind Anti-IL1A Neutralizing antibody (BME100097) in a linear range of 3.05-740.74 ng/mL. In order to specifically detect BME100097, mouse anti-human Fab-specific antibody was used as detection antibody.




Figure 2. Flow cytometry analysis under cell membrane permeable condition with 1 μ g/mL Anti-IL1A (bermekimab biosimilar) mAb (BME100097) on HEK293 cells transfected with Human IL1A protein (Blue histogram) or HEK293 transfected with irrelevant protein (Red histogram).

