

PRODUCT INFORMATION

P01584 Uniprot ID

S-78989;S78989;S 78989;VPM 087;VPM-087;VPM-87;XMA-005.2;XMA-0052;XMA52;XOMA-052;XOMA052;XOMA52;XOMA-52 **Common Name**

Conjugate Unconjugated

IL-1;IL1F2;IL1beta;IL1-BETA Synonyms

Applications N/A Recommended Dilutions N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species N/A lgG2 IgG type Reactivity Human Target IL1B

Anti-IL1B(gevokizumab biosimilar∏IgG2) mAb Description

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

> Email: info@dimabio.com Website: www.dimabio.com

Background Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Usage



