

PRODUCT INFORMATION

Uniprot ID Q9NPF7

LY 3074828;LY-3074828;LY3074828 **Common Name**

Conjugate Unconjugated

Synonyms P19;SGRF;IL-23;IL-23A;IL23P19

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

Host Species N/A IgG4 IgG type Reactivity Human **Target** IL23A

Description Anti-IL23A(mirikizumab biosimilar IgG4) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in

lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store Storage & Shipping

at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

témperature.

Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only

