

PRODUCT INFORMATION

Uniprot ID P01589

BIIB-019;BIIB019;Dacliximab;RO **Common Name**

247375;RO247375;RO-247375;Zenapax;Zinbryta

Conjugate Unconjugated

p55;CD25;IL2R;IMD41;TCGFR;IDDM10;IL2-RA;IL-2-RA **Synonyms**

Applications N/A

Recommended

N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation &

Reconstitution

Host Species N/A IgG type lgG1 Reactivity Human IL2RA **Target**

Description Anti-IL2RA(daclizumab biosimilar IgG1) mAb

Delivery 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid

Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping repeated freezing and thawing). Lyophilized proteins

are shipped at ambient temperature. Research grade biosimilar. Not for use in

Background therapeutic or diagnostic procedures for humans or

animals.

Usage Research use only

