

PRODUCT INFORMATION

|                              |   |
|------------------------------|---|
| Uniprot ID                   | P01589  |
| Common Name                  | AU-007;AU007;AU 007;BD-8;BD8;BD 8;BDG8;BD-G8;BD G8  |
| Conjugate                    | Unconjugated  |
| Synonyms                     | p55;CD25;IL2R;IMD41;TCGFR;IDDM10;IL2-RA;IL-2-RA   |
| Applications                 | N/A   |
| Recommended Dilutions        | N/A   |
| Formulation & Reconstitution | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.  |
| Host Species                 | N/A   |
| IgG type                     | IgG1  |
| Reactivity                   | Human   |
| Target                       | IL2RA   |
| Description                  | Anti-IL2RA(imneskibart biosimilar) IgG1) mAb  |
| Delivery                     | 3-4 weeks   |
| Storage & Shipping           | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Background                   | Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.   |
| Usage                        | Research use only   |

