Cat. No. BME100475



## **PRODUCT INFORMATION**

Uniprot ID P05113

Common Name AQ82742999;AQ-82742999;GSK

294;GSK-294;GSK294;GSK-3511294;GSK3511294

Conjugate Unconjugated

**Synonyms** EDF;TRF;IL-5;Interleukin-5

Applications N/A
Recommended N/A
Dilutions

iutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis

for specific instructions of reconstitution.

Host Species N/A
IgG type IgG1
Reactivity Human
Target IL5

**Description** Anti-IL5(depemokimab biosimilar ☐ IgG1) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

Email: info@dimabio.com Website: www.dimabio.com

**Storage & Shipping** intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only

