Cat. No. BME100644



## PRODUCT INFORMATION

P05113 **Uniprot ID** 

CDP-835;CDP835;CDP 835;CEP-38072;CEP38072;CEP 38072;CINQAERO;CINQAIR;Cinquil;CTx-55700;CTx55700;Ctx 55700;DCP-835;DCP835;DCP 835;SCH 55700;SCH55600;SCH-55700;TRFK-5;TRFK 5;TRFK5 **Common Name** 

Conjugate Unconjugated

EDF;TRF;IL-5;Interleukin-5 **Synonyms** 

**Applications** N/A Recommended **Dilutions** 

N/A

Formulation & Reconstitution

Storage & Shipping

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of

reconstitution.

**Host Species** IgG type IgG4 Reactivity Human **Target** IL5

**Description** Anti-IL5(reslizumab biosimilar IgG4) mAb

**Delivery** 4-6 weeks

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and

Email: info@dimabio.com Website: www.dimabio.com

thawing). Lyophilized proteins are shipped at ambient

temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background** 

Usage Research use only

