

PRODUCT INFORMATION

Uniprot ID P05231

PF-04236921;PF-4236921;PF 04236921;PF 4236921;PF04236921;PF4236921;TOUR-006;TOUR Common Name

Conjugate Unconjugated

CDF;HGF;HSF;BSF2;IL-6;BSF-2;IFNB2;IFN-beta-2 **Synonyms**

Applications N/A

Recommended

N/A **Dilutions**

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of Formulation & Reconstitution reconstitution.

Host Species IgG type IgG2 Reactivity Human **Target**

Description Anti-IL6(pacibekitug biosimilar ☐IgG2) mAb

Delivery

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

> Email: info@dimabio.com Website: www.dimabio.com

Storage & Shipping

temperature.

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. **Background**

Usage Research use only

