Cat. No. BME100552



PRODUCT INFORMATION

Uniprot ID P08887

Common Name BCD 089;BCD-089;BCD089;BCD89;BCD-89;BCD 89;Ilsira

Conjugate Unconjugated

Synonyms IL6Q; gp80; CD126; HIES5; IL-6R; IL6RA; IL6RQ; IL-1Ra; IL-6RA; IL6QTL; IL-6R-1Ra; IL-6RA; IL6QTL; IL-6R-1Ra; IL-6RA; IL6QTL; IL6QTL; IL-6RA; IL6QTL; IL6QTL;

Applications Recommended Dilutions N/A

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. Formulation & Reconstitution

Host Species lgG1 IgG type Human Reactivity IL6R **Target**

Description Anti-IL6R(levilimab biosimilar IgG1) mAb

Delivery 4-6 weeks

Email: info@dimabio.com Website: www.dimabio.com

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. Storage & Shipping

Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals. Background

Usage Research use only

