Cat. No. BME100669



## **PRODUCT INFORMATION**

**Uniprot ID** P05231

CNTO-136;CNTO136;CNTO **Common Name** 136;PLIVENSIA;Shirukumabu

Conjugate Unconjugated

CDF;HGF;HSF;BSF2;IL-6;BSF-2;IFNB2;IFN-beta-2 **Synonyms** 

**Applications** N/A Recommended N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis Formulation & Reconstitution

for specific instructions of reconstitution.

**Host Species** IgG1 IgG type Reactivity Human IL<sub>6</sub> **Target** 

**Description** Anti-IL6(sirukumab biosimilar IgG1) mAb

4-6 weeks **Delivery** 

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

Email: info@dimabio.com Website: www.dimabio.com

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient

témperature.

Research grade biosimilar. Not for use in

**Background** therapeutic or diagnostic procedures for humans

or animals.

Research use only **Usage** 

