

PRODUCT INFORMATION

Uniprot ID P06213

RZ-358;X 358;XOMA

Common Name

247;XOMA-247;XOMA247;XOMA 358;XOMA-358;XOMA358;XPA.15.247

Conjugate Unconjugated HHF5;CD220 **Synonyms**

Applications N/A Recommended N/A

Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before Formulation & Reconstitution lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Host Species N/A IgG type IgG2 Reactivity Human **INSR Target**

Anti-INSR(ersodetug biosimilar IgG2) mAb **Description**

4-6 weeks **Delivery**

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

temperature.

Research grade biosimilar. Not for use in **Background**

therapeutic or diagnostic procedures for humans

or animals.

Usage Research use only



Email: info@dimabio.com Website: www.dimabio.com

