

## **PRODUCT INFORMATION**

**Uniprot ID** Q9UBP8 TBC **Common Name** 

Conjugate Unconjugated

kidney associated antigen 1;kidney-associated antigen 1;MGC78738;RU2 antisense gene **Synonyms** 

protein; RÚ2AS

**Applications** N/A Recommended N/A Dilutions

Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before Formulation & Reconstitution lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Host Species** N/A IgG type IgG1 Reactivity Human KAAG1 **Target** 

**Description** Anti-KAAG1(olintatug biosimilar IgG1) mAb

4-6 weeks **Delivery** 

Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not

intended for use within a month, aliquot and store Storage & Shipping at -80°C (Avoid repeated freezing and thawing).

Lyophilized proteins are shipped at ambient

temperature.

Research grade biosimilar. Not for use in **Background** 

therapeutic or diagnostic procedures for humans

or animals.

**Usage** Research use only





